



Bure Rise & Go™

StandTall walker with power rise



Bure Rise & Go™

StandTall walker with power rise

The Bure Rise & Go is a cost-effective walker enhanced with a patented power rise function. Thanks to the electric power rise function, Bure Rise & Go becomes a combined stand-up, mobility and walking aid all in one single product.



Step 1

The harness is fitted over the patient's lower back and bottom, the straps adjusted and the leg supports placed in the correct position.



Step 2

The patient holds on to the handles and the electrically driven power rise function begins. The patient follows a natural pattern of movement with the aid of the unique harness and the vertical lift action.



Step 3

On completion of the power rise, the leg supports are swung away and the patient is free to walk using the StandTall walker. Alternatively, patients place their feet on footplate and are ready to be moved once power rise is complete.

Unique power rise assistance

Bure Rise & Go saves time, space and money. It acts as a stand-up, mobility and walking aid in one and the same product, which makes the concept extremely cost effective.

The Bure Rise & Go enables better work ergonomics for carers, trainers and anyone else working with patients. Thanks to the smart power rise assistance, standing up becomes a safe, confident and natural action that does not wear out the backs of health care personnel.

The product consists of a walker and its associated harness and straps. Harnesses are available in different sizes (S, M, and L) and they are split at the rear to make e.g. WC visits easier.

When the patient stands up the electric power rise function is employed and the attachment reduces the distance between the patient and the walker until the walker is in the standing position. The Bure Rise & Go contributes to reduced workloads and is an investment in personnel health and the patient's rehabilitation.

Stable

The walker has anatomical cushions that provide relief for shoulders, arms and the neck and shift the focus to larger muscle groups such as leg and chest muscles instead. This improves stability and makes the walker easier to steer.

Flexible stepless adjustment

Electrically driven height adjustment makes it easy to pre-set the walker to the desired height. Once the patient is standing in position the height can be steplessly adjusted for comfort and the proper support for an optimal walking configuration.

Power Rise

Bure Rise & Go is a unique aid for assisting patients or users during the standing manoeuvre. Personnel burdens, and thus repetitive strain injuries, are minimized.

Waterproofing

All Bure Rise & Go electrical components have IP 66 international protection ratings and can therefore be used in wet spaces.

Stirrups for stability

To provide extra safety Bure Rise & Go is fitted with a stirrup to ensure the walker remains stationary during the standing manoeuvre.

Open, reinforced design

The open design provides the same stability but with greater room for manoeuvre. The patient does not feel shut in and it makes assistance with standing up and dressing easier. Bure Rise & Go is also reinforced and extended for extra stability during the standing manoeuvre.

Easily reachable hand controls

The hand controls are located such that both of the carer's or user's hands can easily reach them.

Ergonomic handles

Ergonomically designed handles that are steplessly adjustable in every direction for optimum individual adaptation. This is extremely important for stroke patients who can thus reach their "usual" hand position.

Folding armrests

Bure Rise & Go's anatomically designed armrests can be folded out to provide the user with an easy means of reaching the handles during the standing manoeuvre.

Sustainable

The walker is designed to last a long while. For example it has extra reinforcement bands around the arm rest cushions to prevent them from splitting from impacts to the side.

Unique power rise harness

The power rise harness is designed to reach as far down across the bottom as possible. This provides optimum support which results in a natural standing action.

Lockable castors

Bure Rise & Go has four individually lockable castors as standard.

Leg support for the standing manoeuvre

Bure Rise & Go is fitted with leg supports that make for a better standing manoeuvre; they are height adjustable and can be moved aside simply to create greater space.

Handbrakes and directionally lockable castors

There are one or two optional hand brakes that can be positioned so that the user can reach them from the walking position. The handbrakes act on the rear castors and the directional locks make it easier for the user to stay on course while walking.

Technical specifications

Dimensions apply to walkers with 125 mm castors.

CE Tested and approved according to SS-EN 12182:1999

Part No.		Width, external	Max user weight	Length	Height	Armrests/inside dimension
56-316	Bure Rise & Go StandTall walker	77 cm	150 kg	92 cm	97-137 cm	34-67 cm
56-383	Footplate Rise & Go	-	150 kg	-	-	-
56-384-S/M/L	Harness (available in size S, M, L)	-	-	-	-	-

Bure Rise & Go™

StandTall walker with power rise

- **BURE RISE & GO** is a complete solution that saves time, space and money which means great cost efficiency.
- **SAFE, CONFIDENT** standing and mobility mean a better working environment for health care personnel.
- **THE UNIQUE** power rise assistance provides a natural standing action that allows the patient to participate.
- **BURE RISE & GO** is a patented, protected proprietary GATE Rehab Development AB design.
- **THE WALKER CAN** be used with a footplate that makes patient transits easy.
- **THE HARNESS IS WASHABLE** at 60°, and thanks to its split design it is simple and flexible to wear even during visits to e.g. the WC.
- **BURE RISE & GO** is equipped with a smooth, steplessly variable electric height adjustment.

The footplate
is ideal for
use in patient
transits!




Watch the Bure Rise & Go video at www.gaterd.se



REHAB DEVELOPMENT AB™

Would you like to know more?

Contact us at:
Tel +46 (0)3713 1800
Gate Rehab Development AB
Isberga Säteri,
SE-333 91 Smålandsstenar
info@gaterd.se
www.gaterd.se



GATE COLLECTION